

Ready AAMI COMPLIANT Sets THREE 25-POUND SETS Go IMMEDIATE-USE



The New AMSCO® 400 Series Small Steam Sterilizer.

Compliant. Larger capacity. Immediate-use.

Ready to keep pace with life in any surgical environment.

Compliant. Larger capacity. Ready for immediate-use.

Compliant

- AAMI-ST8:2008
- AAMI-ST79:2010

Processes up to three 25-lb. sets

- Increased shelf clearance allows processing of larger ortho and bariatric sets and large retractors
- Requires no additional floor space

Model	Size	Type	Load
400 Series 16	16 x 16 x 26"	Prevac	2 trays/4 fabric
400 Series 20	20 x 20 x 38"	Prevac	3 trays/6 fabric
400 Series 16S	16 x 16 x 26"	SFPP	2 trays/4 fabric
400 Series 20S	20 x 20 x 38"	SFPP	3 trays/6 fabric

Customizable cycles/cycle names Variety of pre-programmed cycles

Available Cycles:	Sterilize Temp.	Sterilize Time	Dry Time
Immediate-use Prevac	270°F (132°C)	4 Min.	1 Min.
Immediate-use Gravity	270°F (132°C)	10 Min.	1 Min.
Immediate-use Gravity	270°F (132°C)	3 Min.	1 Min.
Prevac	270°F (132°C)	4 Min.	30 Min.
Prevac	275°F (135°C)	3 Min.	30 Min.
Prevac	270°F (132°C)	4 Min.	5 Min.
Prevac	270°F (132°C)	4 Min.	20 Min.
SFPP *	270°F (132°C)	4 Min.	5 Min.
SFPP *	270°F (132°C)	4 Min.	30 Min.
SFPP *	275°F (135°C)	3 Min.	30 Min.
SFPP *	270°F (132°C)	4 Min.	20 Min.
Gravity	250°F (121°C)	30 Min.	15 Min.
Gravity	270°F (132°C)	15 Min.	30 Min.
Gravity	250°F (121°C)	30 Min.	30 Min.
Gravity	270°F (132°C)	25 Min.	15 Min.
Liquid	250°F (121°C)	45 Min.	N/A
Leak Test	N/A	N/A	N/A
DART Test (Bowie Dick)	270°F (132°C)	3.5 Min.	1 Min.
DART Warm-up	270°F (132°C)	3 Min.	1 Min.

* 16S and 20S units only

Safe, intuitive operation

- Hands-free, foot-activated door
- One-touch control with security access code
- Simple programming; customizable cycle names
- Control/countdown screen visible from 15'
- Display messages in complete phrases
- Ink-on-paper impact printer
- Controls at eye level for easy access
- Programmable utilities startup/shutdown
- Water-saving temperature control

Adaptable configurations

- Single Door
- Double Door
- Cabinet Enclosed/Freestanding
- Recessed
- House Steam
- Integral Steam Generator with Flush and Drain
- 208 Volts
- 240 Volts
- 480 Volts

Reliably engineered

- Non-lubricated, steam-activated door gasket with two-year warranty
- Teflon-sleeved shelves with up to 18" clearance
- Easy access modular piping
- Fully jacketed pressure vessel
- 316 L stainless steel construction
- 15-year chamber warranty

Accessories

- Rack/Shelves (16" and 20" units)
- Loading Car/Transfer Carriage (20" only)
- Seismic Tie-Down Kit



STERIS Corporation
5960 Heisley Road, Mentor, OH
44060-1834 • USA
440-354-2600 • 800-548-4873
www.steris.com

